



SAVANNAH

m-oem Raman Module


DROP-IN RAMAN MODULE

- Research & Development
- Clinical Diagnostics
- Process Analytics
- Portable Sensing
- Raw Material ID
- Food and Agriculture
- Semiconductor Characterization



SAVANNAH: Get your technology to market faster with unrivaled Swiss-quality engineering!

Experience the power of innovation with our cutting-edge compact Raman module featuring an integrated laser and spectrograph.



KEY BENEFITS

- Streamline research and development efforts
- Slash development costs
- Speed up your go-to-market timeline
- Proven and widely established Raman platform

KEY FEATURES

- Compact size: 12.32 x 11.22 x 4.71 cm
- Integrated and highly stabilized –785 nm laser w/interlock safety options
- Versatile Sampling Options for liquids, surfaces, powders and slurries.
- Orbital Raster Scanning (ORS)
- Software Development Kit (SDK) included
- Bluetooth and USB communications

